



## Photon Control and Data Systems

# A Simple Online Analysis

## Example

Added by [Christopher O'Grady](#), last edited by [Christopher O'Grady](#) on Sep 18, 2009 ()

Labels:

- **ssh pslogin**
- **ssh atca211**
- **cp -r -d ~cpo/ana .** (this copies the analysis code to your area, only has to be done once)
- **cd ana** (change directory to "ana")
- choose a data file in /reg/d/pcds/amo/onlineo (they are labelled with the time at which they were recorded).

The "user modifiable code" lives in the file myana.cc. This only uses simple arrays, loops, and if statements. To compile, run, and examine the output use the following three commands:

- **./comp** (compiles the code, takes about 2 seconds)
- **./myana -f /reg/d/pcds/amo/onlineo/20090912-060451-o.xtc -n 100** (analyzes 100 events, leave off -n flag to analyze all events)
- **./root** (brings up ROOT gui to examine plots)

The output file is given the same filename as the input, but with a .root extension. The current myana.cc creates a histogram that averages together all shots (called a "profile" histogram) and saves the first 5 shots as separate histograms. There are some comments in myana.cc.

### Common

Electronic Logbook  
XES Common  
LUSI Common  
XTOD  
Supported Control & DAQ Devices

### Instruments

AMO  
CXI  
XPP  
XCS  
LUSI Diagnostics  
SXR  
MEC  
CAMP  
EBIT

### Detectors

CXI 2D Detector  
XPP-2D-Detector  
pnCCD-Detector

### Subsystems

Networking  
Data Acquisition  
Data Management  
Controls  
Laser Control  
Accelerator Timing  
System

Femto-Second Timing  
Test-Stands

Personel Protection  
System

Machine Protection  
System

XES Laser Safety  
Systems

User Safeguards

## Tools

Confluence How-To

Jira How-To

Internal

## Support

Controls How-To

[Controls Problems &  
Solutions]

## Misc

Misc



Page Operations



Browse Space

## Within the root GUI:

- double-click on the "ana" folder. this is the directory where the analysis output lives.
- double-click on the .root file that was output from the analysis (same filename as input, but with .root extension)
- that one root file is then selected. unfortunately, you have to double-click on it again.
- double-click on the histogram you wish to view. You can zoom in on axes by clicking and dragging on the axis, for example.
- you can quit root by typing ".q" at the terminal prompt

More information on root can be found at  
<http://root.cern.ch>

Powered by [Atlassian Confluence](#) 2.10.2, the [Enterprise Wiki](#). [Bug/feature request](#) – [Atlassian news](#) – [Contact administrators](#)